

# Goggles Frame

# FHS/T00



**MG-1.8; MG-2.3; MG-2.7**

The magnifier lenses are compatible with both the **FHS/T00** goggles and the **FHS/T01** miner's lamp, providing 1.8x, 2.3x and 2.7x magnification.


These goggles provide emission filters for fluorescent protein visualization conveniently located on the observers' head. They are designed to use the **FBL/Basic-B&N-01 Excitation Stand** as a light source, but can be used with the **CVS 8 Camera/Video Stand** and the **UFP-MDS-cc/HB Flashlight**, too.

The emission filters are available for the typical FPs: ECFP, EGFP, EYFP, RFP and OFP. For custom emission filters please contact us.

The emission filters are interchangeable, so that the same goggles can be used with different kinds of FPs.

The MG-1.8, MG-2.3 and MG-2.7 magnifier lenses can be attached to the emission filters.

The goggles and the emission filters are sold as separate units.



Emission Filter Model Number	Fluorescent Proteins	Open (Wavelength nm)
FS/TEF-1C2	Cyan	480-505 nm
FS/TEF-2G2	Green	500-515 nm
FS/TEF-3GY1	Green & Yellow	500-650 nm
FS/TEF-3GY2	Green & Yellow	500-550 nm
FS/TEF-4Y1	Yellow	520-650 nm
FS/TEF-4Y2	Yellow	520-550 nm
FS/TEF-4Y3	Yellow	520-540 nm
FS/TEF-4R2	Red	590-660 nm
FS/TEF-4O2	OrangeRed	557-590 nm



BIOLOGICAL LABORATORY EQUIPMENT  
MAINTENANCE AND SERVICE LTD.

[www.bls-ltd.com](http://www.bls-ltd.com)

Újszász u. 47/E, H-1165 Budapest, HUNGARY

Tel.: +36-20-4297735 E-mail: [bls@bls-ltd.com](mailto:bls@bls-ltd.com) Website: [www.bls-ltd.com](http://www.bls-ltd.com)