

Guide For FHS/T00 and FHS/T01



FHS/T00



FHS/T01



MG - 18; MG - 23; MG - 27





MG - 18; MG - 23; MG - 27

MAGNIFIER LENSES

EMISSION FILTER - FS/TEF-xxx

| Emission Filter Model Number | Fluorescent Proteins | Open (Wavelength nm) |
|------------------------------|----------------------|----------------------|
| FS/TEF-1C2 | Cyan | 480-505 nm |
| FS/TEF-2G2 | Green | 500-515 nm |
| FS/TEF-3GY1 | Green & Yellow | 500-650 nm |
| FS/TEF-3GY2 | Green & Yellow | 500-550 nm |
| FS/TEF-4Y1 | Yellow | 520-650 nm |
| FS/TEF-4Y2 | Yellow | 520-550 nm |
| FS/TEF-4Y3 | Yellow | 520-540 nm |
| FS/TEF-4R2 | Red | 590-660 nm |
| FS/TEF-4O2 | OrangeRed | 557-590 nm |



**Camera/Video Stand
CVS 8**



**Flashlight
UFP-MDS-CC/HB**

**Excitation stand
FBL/Basic-B&N-01**

LIGHT SOURCE FS/ULS-02-cc



| Light Source Model Number | Excitation Wavelength | Emission for Fluorescent Proteins |
|---------------------------|-----------------------|-----------------------------------|
| FS/ULS-02RB | 440-460 nm | Cyan |
| FS/ULS-02B2 | 460-495 nm | Green & Yellow |
| FS/ULS-02G2 | 515-555 nm | Red |
| FS/ULS-02OR | 567-588 nm | OrangeRed |



BIOLOGICAL LABORATORY EQUIPMENT
MAINTENANCE AND SERVICE LTD.

www.bls-ltd.com

Újszász u. 47/E, H-1165 Budapest, HUNGARY

Tel.: +36-20-4297735 E-mail: bls@bls-ltd.com Website: www.bls-ltd.com